

## Chromatin IP

Use 3-10cm plates per IP.

### Crosslinking:

Crosslink by adding 270 $\mu$ L 37% formaldehyde to media in 10cm plate to achieve 1% FMA. Incubate cells at room temp, 10min while gently shaking.

Wash cells 2X w/ 5mL ice-cold PBS.

Add 0.5mL rinse buffer to each plate and scrape cells off plate and place in eppy. Incubate cells 15min @ 30°C, and spin down @ 2000g, 5min, RT.

- Wash and combine pellets w/ 1mL ice cold PBS, and spin down @ 2000g, 5min, 4°C.
- Lyse cells w/ 1mL buffer I. Incubate on ice 5min, and spin down @ 2000g, 5min, 4°C to collect nuclei.
- Wash nuclei w/ 1mL buffer II, and spin down @ 2000g, 5min, 4°C.

Lyse nuclei w/ 2mL SDS lysis buffer. Can store nuclear extract @ -80°C or proceed to sonication.

### Sonication:

Sonicate using T. Johnson sonicator @ 30% power w/ 10s pulses in a FACS tube. Allow nuclear extract to rest 1min on ice between pulses. Sonicate a total of 3min.

After sonication, spin nuclear extract 10min @ 14,000g, 4°C.

Dilute samples 1:10 in dilution buffer, aliquot 100 $\mu$ L for input. Can store diluted nuclear extracts @ -80°C or proceed to IP.

### Immunoprecipitation:

Use 5mL extract for each IP. Preclear extract 2X 1hr w/ 100 $\mu$ L blocked CL-4B beads while rotating. Remove supernatant from beads and incubate overnight w/ antibody (usually 1-2 $\mu$ g antibody).

### Blocking Beads:

Wash slurry at least 3X w/ 1mL dilution buffer to remove EtOH. Add 20 $\mu$ g sonicated salmon sperm DNA per mL of slurry. Rotate at 4°C for 1h. Store beads at 4°C.

The following morning add 40 $\mu$ L of block protein G slurry and rotate 2h @ 4°C.

### Washes and Elution:

Spin down IP and remove all supernatant. Use first TSE I wash to transfer beads into eppy.

- Wash 2X 15min w/ 1mL TSE I while rotating at RT.
- Wash 1X 15min w/ 1mL TSE II while rotating at RT.
- Wash 1X 15min w/ 1mL buffer III while rotating at RT.

Remove all traces of wash buffer by rinsing beads 3X w/ 1mL TE.

Elute 3X w/ 100 $\mu$ L of freshly prepared elution buffer. Rotate in elution buffer 5min @ RT, then vortex 30s, and spin down to collect eluate. Combine eluate (300 $\mu$ L total) and then spin down again to remove all beads.

### Reverse Crosslinks:

Add 200 $\mu$ L elution buffer to 100 $\mu$ L input. Incubate input and eluate @ 65°C for 4-6h to reverse crosslinks.

After crosslinks are reversed, add 30 $\mu$ L 3M NaOAc pH 5 and 660 $\mu$ L 100% EtOH, and precipitate o/n @ -20°C.

Next day, spin down DNA 30min @ 14,000g, 4°C. Resuspend in 38 $\mu$ L TE, add 10 $\mu$ L PK buffer and 2 $\mu$ L (100 $\mu$ g) proteinase K. Digest @ 42°C, 1-2h. Bring up volume to 300 $\mu$ L w/ TE, and add equal volume of phenol chloroform. Vortex 15s and spin @ 14,000g for 5min @ RT.

Collect 280 $\mu$ L of aqueous layer and precipitate overnight w/ 28 $\mu$ L 3M NaOAc pH 5.0, 626 $\mu$ L 100% EtOH, and 10 $\mu$ g glycogen.

Next day, spin down DNA. Resuspend DNA in 10mM Tris pH 8.5. Use 200 $\mu$ L for input and 20 $\mu$ L for ChIP samples.

Perform PCR w/ 1 and 5 $\mu$ L of DNA.

Rinse Buffer: 100mM Tris pH 8.8 10mM DTT	Buffer II: 200mM NaCl 10mM EDTA 10mM Hepes pH 6.5	Dilution Buffer 1% triton 2mM EDTA 150mM NaCl 20mM Tris pH 7.4 1mM PMSF	TSE II: 0.1% SDS 1% triton 2mM EDTA 20mM Tris pH 7.4 500mM NaCl
Buffer I: 0.25% triton 10mM EDTA 0.5mM EGTA 10mM Hepes pH 6.5	SDS lysis buffer 1% SDS 10mM EDTA 50mM Tris pH 7.4 protease inhibitors	TSE I: 0.1% SDS 1% triton 2mM EDTA 20mM Tris pH 7.4 150mM NaCl	Buffer III: 0.25M LiCl 1% NP40 1% Deoxycholate 1mM EDTA 10mM Tris pH 7.4
TE: 10mM Tris pH 7.4 1mM EDTA	Elution buffer: 1% SDS 0.1M NaHCO <sub>3</sub>		